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
Operating Instructions

1 PURPOSE AND DESCRIPTION

The Dipole system is tuned to the Phywe microwave transmitter. It consists of a base plate which is a non-conductor of electricity and has several quadratic arrangements of conical holes in which metal rods can be fitted.

When these metal rods are irradiated by the microwave transmitter they act as secondary excited transmitters, whose interfering wave field can be examined for various geometrical arrangements of the rods.

Going out from the central hole in the base plate, the holes in the middle of the sides of the quadratic arrangements are each at a distance of $\lambda/4$ or $\lambda/2$ from each other (the length of the transmitted primary wave is $\lambda = 3.2$ cm).

2 HANDLING

The 9 metal rods, whose length is about four times the wavelength, each have a conically shaped end which fits firmly into the holes in the base plate. The other end is coloured white for easy recognition of the type of arrangement of the inserted rods. The length of the base plate rod has been so chosen that the middles of the dipoles are at the same height as the microwave transmitter and receiver when all equipment is held, for example, by tripods on the experimental surface.

When not in use, the rods can be stored in clamps underneath the base plate.

3 EXPERIMENTS

Experiment 1: Fit 5 rods on the base plate so that one rod is in the central hole, and the other four are one in each of the four middle holes of the fourth quadrate out from the central hole. The distance between the rods in one row is

now one wavelength. Now transmit a sound frequency modulated wave in the direction of a rod row, with the microwave transmitter situated approximately 20 cm. away. Position the microwave receiver, connected to an amplifier and a loudspeaker, perpendicularly to the irradiated rod row. The receipt of a signal can be acoustically detected, but on the slightest swivelling of the base plate, the signal disappears. The rod system therefore radiates at right angles to the exciting wave. The rods act as transmitter, swinging in phase because of their spatial arrangement. This phase equality in time is also spatially only maintained in preferred directions.

Experiment 2: Add 4 further rods to those used above, one in each of middles of the sides of the second quadrate. The distance between the rods in each row is now half a wavelength; the signal perpendicular to the direction of the primary wave disappears almost completely.

The rods are now excited alternately by opposite phases, because of the difference in transmission time of the primary waves. This phase situation is maintained at the position of the receiver, because the path lengths are hardly different.

Experiment 3: Insert all 9 rods in a single line. The primary wave meets this line of rods at any incident angle. The momentary phase differences in the currents generated in the individual rods are given by the differences in the paths, which are dependant on the incident direction of the primary wave. Because of the path differences of the emitted secondary waves, phase equality, and so a receivable signal, is only given in the direction of the emitted angle. The line of rods acts in its reflection behaviour like a metal plate.

4 LIST OF EQUIPMENT

Dipole system
Microwave power supply, 220 VAC
Microwave transmitter with klystron
Microwave receiver
LF amplifier, 220 V
Loudspeaker, 8 Ω/5 kΩ
Tripod base -PASS-

06860.01
11740.93
11740.01
11740.02
13625.93
13765.00
02002.55

5 NOTE ON THE GUARANTEE

We guarantee the instrument supplied by us for a period of 6 months. This guarantee does not cover natural wear nor damage resulting from improper handling. The manufacturer can only be held responsible for the function and safety characteristics of the instrument, when maintenance, repairs and changes to the instrument are only carried out by the manufacturer or by personnel who have been explicitly authorized by him to do so.

Microwave power supply

11740.93

Operating instructions

The unit complies with the corresponding EC guidelines.

1. PURPOSE AND DESCRIPTION

The microwave power supply serves to operate the microwave transmitter with klystron. Within the bounds of a transmitting and receiving system, the properties of electromagnetic wave fields, such as propagation, reflection, refraction, interference, diffraction and polarization can be examined, using centimetre waves as an example. Further to this, technical uses such as radar and directional radio can be demonstrated.

The microwaves can be radiated either without modulation or amplitude modulated with low frequency voltage. The radiation diagram has directional characteristics.

2. PARTS AND OPERATING ELEMENTS

2.1 The power supply

The following operating and functional elements are to be found on the front plate (Fig. 1):

1. Mains switch with control light

2. Set of buttons for modulation

With button 2.1, the form of the modulation signal is determined. With button 2.2 the choice of either internal or external modulation voltage is made. The modulation types which are possible are given in the Table:

Button 2.1	Button 2.2	Modulation by
0	0	ext. signal voltage voltage peaks cut off
1	0	ext. signal voltage in unchanged form
1	1	int. sinus voltage 50 Hz
0	1	int. rectangular voltage 50 Hz

(button not pressed = 0; button pressed = 1)

When there is no signal voltage on the pair of sockets 6, then, with button 2.2 not depressed, unmodulated microwaves are radiated.

3. Adjusting knob for reflector voltage

For the optimum adjustment of the reflector direct voltage of the reflex klystron, to achieve the maximum radiation intensity of the microwaves.

When there is no signal voltage on the pair of sockets 6, then, with button 2.2 not depressed, the microwaves are not modulated.

4. Klystron socket

Multi-pole socket for connection of the microwave transmitter.

5. 4 mm socket

with casing mass, connected to the earthed socket of the electrical supply.

6. Input of external modulation

Pair of 4 mm sockets for the feed-in of external signal voltages ($U_{\sim} = \text{max. } 4 \text{ V}$; $f = 50 \text{ Hz} \dots 10 \text{ kHz}$).

2.2 Microwave transmitter with klystron

The microwave transmitter is equipped with a reflex klystron as oscillation generator. The transmitter radiates the waves out for bundling over a horn radiator of rectangular cross-section. The direction of oscillation of the electrical vector of the polarized microwave runs parallel to the narrow side of the funnel.

For a defined variation of the direction of oscillation in experiments, the stem of the transmitter can be screwed in in three directions which are perpendicular to each other.


Fig. 2: The transmitter in diagram form

- 1 Transmitting tube (klystron)
- 2 Hollow space resonator
- 3 High frequency lead
- 4 Reflection seal of the wave guide
- 5 End of the lead for exciting the wave guide
- 6 Wave guide
- 7 Direction of oscillation of the electrical field vector
- 8 Direction of emergence of the wave field

3. THE MICROWAVE RECEIVER

For the receiving of the microwaves there is a choice between the microwave receiving dipole and the microwave directional receiver, both of which use a silicon diode as receiving antenna. The voltage generated by the microwaves is measured at this, either directly or over an amplifier. The receivers have a BNC socket for their connection.

Wave fields can be measured out almost as points with the microwave receiving diode, whereas with the directional receiver, the radiation is collected over a larger angular area and so delivers a higher voltage than the dipole when they are at the same distance from the transmitter. Two threaded borings enable the receiving antenna to be held in either of two positions, which are at 90° to each other.


Fig. 3: The directional receiver in diagram form

- 1 Direction of wave entry
- 2 Wave guide
- 3 Receiving antenna (silicon diode)
- 4 Reflection seal of the wave guide
- 5 to the BNC socket (measurement cable)
- 6 Spring diode connection to the measurement cable
- 7 Insulated passage to the diode
- 8 Direction of oscillation of the electrical field vector

4. MEASURING EQUIPMENT

The receiving antenna delivers a direct voltage signal with non-modulated microwaves, whose size is a measure of the microwave intensity received in the experimental set-up concerned. With amplitude modulation, the voltage given up from the diode changes with the modulation frequency. The following are suitable for the measurement of the diode voltage;

- an oscilloscope
- a mV direct voltage meter
- a low frequency amplifier (only for modulated microwaves)

5. HANDLING

5.1 Starting up for the first time

Plug the klystron into the socket on the microwave transmitter. A socket guide stud and corresponding slot in the socket avoid any possibility of wrong connection. Fit on the connecting cap.

5.2 Assembly of the microwave apparatus

Plug the connecting cable of the microwave transmitter into the socket 4 of the microwave power supply. The transmitter begins to operate about 20 seconds after being switched on, and reaches its stable operating condition after about 3 minutes.

Select the type of radiation with the set of buttons 2 and, if appropriate, put an external modulation voltage on the input 6. Connect the microwave receiver with the measuring equipment by means of a BNC cable. Should this only have 4 mm input sockets, then additionally use a BNC socket / 4 mm plug pair adapter. The receiver should be on the axis of symmetry of the radiated field.

Generally, the transmitter and the directional receiver are used in their normal position (vertical position of the narrow sides of the funnel). The dipole must be aligned vertically to receive with the transmitter in the normal position. To reduce interfering reflections, the transmitter and the receiver should not be held too near to the surface of the table.

Adjust the reflector voltage so that the maximum radiation intensity is given, this can be the case at two different positions of the adjusting knob.

When using an NF amplifier, this should not stand too near to the receiver, as otherwise self-distortion can occur. A loudspeaker or a voltmeter can be connected to the appropriate output of the NF amplifier. As a result of unavoidable scattering and interfering effects, as well as the formation of standing waves, the ideal course of the field can be subject to interference, which is less noticeable with the acoustic signal than with the voltmeter display.

5.3 Changing the silicon diode

The receiver diode is very sensitive to too high voltages. Take care therefore when putting in a new diode, that you yourself are not too highly statically charged. To avoid static charge, touch the connected transmitter with your hand immediately beforehand.

The diode has a metal cap on one of its connecting pins, which must be left on when the diode is used with the directional receiver, but taken off when it is used with the dipole.

6. THE FOLLOWING CONDITIONS MUST BE HELD TO WHEN OPERATING THE TRANSMITTER

1. The transmitter must only be operated at the factory adjusted frequency of 9450 MHz (tolerance +/- 75 MHz). Connection to other components must not cause a change in the frequency.
2. The transmitter must only be used in closed rooms. When a directional antenna is used, radiation towards windows, doors and outer walls must be avoided. The transmitter must not be used in wooden buildings.

7. SPECIFICATIONS

Frequency (9.450 +/- 0.075) GHz
Emitting power 25 mW

Please note:

Microwave transmitters built before 1999 cannot be operated with the newer (built after 1999) microwave power supplies, due to technical modifications made for safety purposes. Our service department offers an economical adaptation of such transmitters.

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9. ACCESSORIES

Microwave transmitter with klystron	11740.01
Microwave receiver	11740.02
Microwave receiving dipole	11740.03
LF amplifier, 220 V	13625.93
Loudspeaker, 8 Ohm / 5kOhm	13765.00

Replacement parts

Klystron	06864.00
Diode 1 N 415 C	06862.01